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Mettler-Toledo  
<http://www.mt.com>

METTLER TOLEDO INSTRUMENTS (SHANGHAI) CO., LTD.  
Printed by MTCS Marcom 2009/01 12107056B

## General information

Chloride electrode is intended for measuring chloride ion concentrations and activities in aqueous, partially aqueous organic solutions. Its preferred applications are to be found in the fields of ground-water, drinking-water, sewage, medicine, pesticide, beverage, food products and soil samples, as well as that of biological fluids (such as urine and perspiration, etc.).

## Specifications

- Type of ion selective membrane:  
NaCl poly-crystal
- Membrane impedance: < 50MΩ
- Measurement range:  
 $5 \times 10^{-5}$  to 1 mol/l Cl<sup>-</sup>
- Operating temperature range: 0 ~ 80° C
- Optimum pH range: 2~11
- Response time:  
when passing from lower to higher concentration sbelow  $10^{-3}$  mol/l Cl<sup>-</sup>, about 1 minute, above  $10^{-3}$  mol/l Cl<sup>-</sup> less than 1 minute; when passing from higher to lower concentrations, several minutes.
- Interfering ions:  
S<sup>2-</sup> must be absent; the molar proportions of the following ions to the chloride ion must lie below the following levels(ratios):  
CN<sup>-</sup>< $2 \times 10^{-7}$ , I<sup>-</sup>< $5 \times 10^{-7}$ , Br<sup>-</sup>< $3 \times 10^{-3}$ , S<sub>2</sub>O<sub>3</sub><sup>2-</sup><0.01, NH<sub>3</sub><0.12, OH<sup>-</sup><80

## Preparatory operations

### a. Activation

Before the chloride electrode was used, it should be soaked in 0.001 mol/l chloride solution for 1 hour and rinsed repeatedly with deionized water till blank potential is steady. The chloride electrode should be dabbed off with paper tissue.

### b. Preparation of standard solution

Dry AR grade sodium chloride at 120° C for 2 hours. Weight out 5.844g NaCl (for a 0.1 mol/l Cl<sup>-</sup> solution), or 1.648g NaCl (for a 1.0 g/l Cl<sup>-</sup> solution), dissolve in deionized water in a 1000ml calibrated flask and make up to the mark.

Prepare solutions with different chloride concentrations by diluting the standard stock solution with deionized water as appropriate.

### c. Reference electrode

Remove cap from electrolyte filling port. Slightly loosen the ground-joint diaphragm (by careful turning) so that out reference electrolyte flows out. Tighten up ground joint again (but not too tight). Top up with 1 mol/l KNO<sub>3</sub> electrolyte to filling port for out-mediate electrolyte and recover the cap. Clutch the electrode and shake up and down in order to reduce air bubbles.

## Calibration and measurement

If you are using **ionic meter**, proceed as follows to measure chloride ion concentration:

### a. 2 point calibration

Prepare two standards; one standard has 10 times difference in concentration from another one. The standard solutions can be prepared with any concentration unit.

According to the menu setting process of your ionic meter, set the lower concentration standard as the first calibration point, and the higher one as the second calibration point.

To 100ml the lower concentration standard add 2ml ISA solution, and stir. Place the chloride electrode and reference electrode together in the mixed solution, start the first point calibration.

To 100ml the higher concentration standard add 2ml ISA solution, and stir. Place two electrodes in the mixed solution, start the second point calibration.

After 2 points calibration, the meter automatically determines the calibration slope.

### b. Measurement Sample

Rinse electrodes with deionized water. Take 100ml sample in a sample container; add 2ml ISA solution and stir. Place two electrodes in the mixed sample solution and start measurement.

If you are using **pH/mV meter**, proceed as follows to measure chloride ion concentration:

### Plotting calibration diagram

Select two or more standard solutions including the sample's chloride concentration. The standard solutions can be prepared with any concentration unit. According to pH/mV meter instruction manual, select mV mode.

To 100ml standard solutions add 2ml ISA solution and stir, respectively. Place the chloride electrode and reference electrode together in the mixed solution and measure the mV of the mixed solution from lower to higher concentrations.

Plot measured mV values on a calibration diagram against the logarithm of the chloride concentration of the corresponding solution (use log graph paper).

### b. Measurement Sample

Rinse electrodes with deionized water after each measurement.

To 100ml sample solutions add 2ml ISA solution and stir. According to above method measure the mV of sample solution. Use the calibration diagram to read off the chloride concentration of the sample solutions from the measured mV values. It is advisable when carrying out lower chloride determination to apply a standard addition method.

## Maintenance

If deposits and scratches were found on the crystal, it can be removed by polishing with a suitable abrasive paste and carefully wiping with alcohol cotton ball in order to renew the sensitive surface of the single crystal. Then the chloride electrode should be reactivated in 0.001 mol/l chloride solution for an hour.


## Storage

After measurement, chloride electrode must be rinsed. The chloride electrode should be stored dry in the air, preferably in the protective tube supplied with the electrode.

## Accessories

LE302 Ref. electrode  
(Order No. 12107202 Or 12107204)  
Bridge electrolyte: 1M KNO<sub>3</sub>  
Ionic Strength Adjuster (ISA): 5M NaNO<sub>3</sub>

Order No.: 12107070 or 12107080

- 1 Electrode cap
- 2 Plastic electrode body
- 3 Internal lead out
- 4 Electric screw
- 5 Ion-selective membrane
- 6 Filling port for inner electrolyte
- 7 Filling port for outer electrolyte
- 8 Internal reference electrolyte
- 9 Ag/AgCl reference element
- 10 Porous ceramic diaphragm
- 11 Outer reference electrolyte
- 12 Ground-joint diaphragm